



REF-LINE®
Look and feel natural.



PoliDent



Products Catalogue



PoliDent



Artificial acrylic teeth

01 REF-LINE

High end anatomically shaped multi-layer tooth line for new generation prostheses

Natural colour and shape

- natural morphology of tooth surface
- incisal cracks and minor irregularities
- transparent incisal area
- visible mamelons in the dentin area
- intensive chromaticity of the gingival region
- including cervical wider shapes of tooth designs for covering metal structures in implants and partial prostheses
- obscured fissural zones on the occlusal surface
- anatomically designed pronounced cusps assure the main taining of the natural tooth appearance and structure in individual matching with the antagonists

Brand new production technology

- perfection of tooth design
- retention opening on the lower part of transcanine teeth provides for easier and faster work

New material




- 50% improvement in tooth wear*
- 15% improved hardness*
- elasticity remains unchanged*

*compared to the existing material




REF-LINE®
Look and feel natural.

Q1 REF-LINE

		
M1	V1,V3	Q1
M3	V3,V4,U3	R2,R4,Q1,Q3,F2,F4
M4	V4,U3,U4	R2,R4,Q3,F2,F4
M6	V4,V6,V8,U4	R4,R6,Q5,F4,F6
N3	V3,V4,U3	R2,R4,Q3,F2,F4
N4	V4,V6,U3,U4	R2,R4,Q3,F2,F4
N6	V6,V8,U4,U6	R4,R6,Q5,F4,F6
W3	V3,V4,U3	R2,R4,Q3,F2,F4
W4	V4,V6,U3,U4	R2,R4,Q3,F2,F4
W6	V6,V8,U4,U6	R4,R6,Q5,F4,F6
H3	V3,V4,U3	R2,R4,Q3,F2
H4	V4,V6,U3,U4	R2,R4,Q3,F2,F4
H6	V6,V8,U4,U6	R4,R6,Q5,F4,F6

ODTENKI SHADES NIJANSE


BL1 - BL3

A1 - A2 - A3 - A3,5 - A4

B1 - B2 - B3 - B4

C1 - C2 - C3 - C4

D2 - D3 - D4



M1

a 40,8
b 8,7
c 7,5

M3

a 43,2
b 9,7
c 8,2

M4



a 45,3
b 9,9
c 8,5

M6

a 49,1
b 10,5
c 9,1



Artificial acrylic teeth



N3

a 43,6
b 9,9
c 8,3



H3

a 43,0
b 10,7
c 7,8



N4

a 45,4
b 10,9
c 8,6



H4

a 44,9
b 11,2
c 8,2



N6

a 48,8
b 11,4
c 9,3



H6

a 48,3
b 12,0
c 8,8



W3

a 43,4
b 9,3
c 8,2



W4

a 45,3
b 9,7
c 8,5



W6

a 48,7
b 10,6
c 9,2





01 REF-LINE



REF-LINE®
Look and feel natural.

01 REF-LINE

Artificial acrylic teeth



V1

a 29,8
b 8,0
c 4,3



V3

a 31,9
b 8,3
c 4,7



V4

a 33,7
b 8,7
c 4,9



V6

a 35,8
b 9,2
c 5,3



V8

a 38,6
b 10,5
c 5,6




U3

a 34,1
b 8,4
c 5,1



U4

a 35,9
b 8,7
c 5,4



U6


a 38,3
b 9,4
c 5,7



REF-LINE®
Look and feel natural.

Artificial acrylic teeth

01 REF-LINE



R2

- a 30,6
- b 32,8
- c 10,6
- d 9,8



R4

- a 33,2
- b 35,4
- c 11,4
- d 10,7



R6


- a 35,0
- b 38,0
- c 12,3
- d 11,3



REF-LINE®
Look and feel natural.

01 REF-LINE

Artificial acrylic teeth



F2

- a 30,4
- b 32,4
- c 10,2
- d 9,5



F4

- a 32,9
- b 35,1
- c 11,1
- d 10,2



F6


- a 35,4
- b 37,7
- c 11,9
- d 11,0



REF-LINE®
Look and feel natural.

Artificial acrylic teeth,

01 REF-LINE



Q1

- a 29,1
- b 31,6
- c 10,2
- d 9,4



Q3

- a 31,9
- b 34,6
- c 11,2
- d 10,3



Q5

- a 33,8
- b 36,7
- c 11,9
- d 10,9



R



Q

extended cervical on premolars



REF-LINE®
Look and feel natural.